






















Product List

	Name	SKU	Description	Contents
	Basic Weekly Testing Service	Y52	Spore Testing services are available for Steam/Autoclave, Chemical Vapor, Dry Heat and Ethylene Oxide sterilization cycles.	52 Tests with Matching Controls, Online Test Results
	Biweekly Testing Service	Y26	Spore Testing services are available for Steam/Autoclave, Chemical Vapor, Dry Heat and Ethylene Oxide sterilization cycles.	26 Tests with Matching Controls, Online Test Results
	Monthly Testing Service	Y12	Spore Testing services are available for Steam/Autoclave, Chemical Vapor, Dry Heat and Ethylene Oxide sterilization cycles.	12 Tests with Matching Controls, Online Test Results
	Weekly Testing Service Kit	Y52P	Spore Testing services are available for Steam/Autoclave, Chemical Vapor, Dry Heat and Ethylene Oxide sterilization cycles.	Kit including: 52 Tests with Matching Controls, Online Test Results, Log Book and 1 package of CI-SSW Steam Integrators (250/package)
	Mesa Record Keeping Booklet	SML MC		1 Booklet
	Smart-Read EZTest	SEZS/5	This unique system uses a real biological indicator (BI) which is incubated, evaluated, and documented in one simple, automated operation.	100 Count 10 Hour Incubation
	Smart-Read EZTest	SEZS/525	This unique system uses a real biological indicator (BI) which is incubated, evaluated, and documented in one simple, automated operation.	25 Count 10 Hour Incubation
	EZTest Steam	EZS/5	EZTest is easy-to-use and produces visual results within 24 hours.	100 Units 24 Hour Incubation
	EZTest Steam	EZS/525	EZTest is easy-to-use and produces visual results within 24 hours.	25 Units 24 Hour Incubation
	EZTest 1410 Incubator	I1410	Compatible with EZTest, adjustable temperatures	

Mesa Canada Product List

	Name	SKU	Description	Contents
	EZTest Crusher	I-CMP	Designed to offer consistent activation of EZTest	1 crusher
	Twindicator	TWN250	The white indicator line on the Twindicator sterilization indicator strip changes to black and provides visible proof that saturated steam has fully penetrated the pack or pouch within which the strip was placed.	Sold in a box of 250 8" strips
	Dry Heat Indicator Labels	CI-DHI	Dry Heat Indicator Labels clearly distinguish processed from unprocessed items when used as directed.	Dry Heat Indicator Labels (Roll of 1000)
	Chemical Indicator Pen	CI-PEN	Ink Changes from Purple to Green when Sterilized!	Class 5 grade ink.
	ProChem SSW Steam Sterilization Integrator	CI-SSW	ProChem SSW Steam Sterilization Integrator is a single use Class 5 Integrating Indicator which complies with ANSI/AAMI/ISO 11140-1:2005 and ANSI/AAMI ST60:1996.	Class 5 Integrator, 250 units per bag
	ProChem SSI	CI-SINT	Use ProChem SSIs to verify adequate steam penetration into every pack/pouch.	Class 4 250 indicators/box
	ProChem EXT	CI-EXT	ProChem EXT is a Class 5 Integrating Indicator (Integrator) that is designed to accurately depict the lethality of steam sterilization.	ProChem EXT is packaged with 250 integrators per box
	Bowie-Dick Green Card Refills	CI-BD134-30	This easy-to-use test requires less storage space than other Bowie-Dick products and is environmentally friendly! Can be used as a permanent record – test information can be recorded on the back.	30 Tests per bag
	Bowie-Dick Green	CI-BD134	This easy-to-use test requires less storage space than other Bowie-Dick products and is environmentally friendly! Can be used as a permanent record – test information can be recorded on the back.	120 per case

Mesa Canada Product List

	Name	SKU	Description	Contents
	Bowie-Dick Green Card Starter kit	CI-BD134SK	The most cost-effective method to perform daily Bowie-Dick testing of pre-vacuum sterilizers.	One holder and 30 test cards
	BD Green Card Holder	CI-BD134TH	Reusable aluminum holder	1 Holder
	Bowie-Dick Record Keeper	LRB-BDSK		50 sheets (10 tests/ sheet)
	Green Card Bowie-Dick Protector Sheets	LRB-BDGC		10 holders per sheet, 50 sheets/box
	Bowie-Dick Record Keeper	LRB-BDSK		50 sheets (10 tests per sheet)